

bravo[®]



WAARDENBURG
MEDICAL





bravo®



The smart and agile companion for small devices.

With its compact and light-weight design, bravo® is the ideal cart for small devices. No matter where it is put to use, this universal companion offers a sturdy basis for a multitude of applications – from the doctor's surgery to the cosmetic studio. The mobile, light-duty equipment cart is also ideal as a carrier system for presentation devices in offices or trade fairs or for testing equipment in industry.

The bravo® equipment cart series stands for dependable brand quality at a reasonable price. It is protected with a European design patent.

Are you looking for a cart that's the perfect fit? On the next pages, we'll show you a range of particularly popular configurations for our bravo®. No more and no less. The bravo® truly is a unique piece that can be put together and equipped completely according to your needs.



Learn about the bravo® in more detail.

Features

- Single-column equipment cart system, light-duty.
- Two-part, single-column aluminium carrier profile back and front, with practical cable ducts.
- Ideal for attaching various shelves, socket strips and diverse, easily adaptable accessories.

Potential function

Computer cart | tablet cart | ergometry cart | hygiene cart | dental cart | projector cart ...

Potential utilisation

Doctors' surgeries | dental surgeries | labs | clinics | beauty studios | consultancy rooms | industrial companies ...

The bravo® system components

- Sturdy cast-aluminium chassis, powder-coated in basalt grey (RAL 7012).
- High-quality double castors as standard (Ø 75 mm) with two-wheel locks and grey wall guards.
- Electrically conductive castors available for a surcharge.
- Available as standard in two different heights (780 mm, 900 mm) or in custom dimensions.
- Two-part aluminium carrier profile, anodised and with integrated cable duct.
- High-quality, sturdy shelf, in part with moulded pull-grip or lateral push-handles.
- Shelves in various sizes, metal or ABS plastic.
- Practical, white, pull-out storage tray (polycarbonate) with storage shelf for dust-free storage of small accessory components.
- Many practical accessories, can be retrofitted at any time, for attaching back and front to the profiles.
- Special editions, custom dimensions and custom colours available on request.
- Balancing weight plate depending on the configuration of the equipment cart system to increase stability optionally available (Ref. No. 27209).

Note

Can be produced for equipment manufacturers customised to requirements with special shelves, holders, components in custom dimensions and custom colours.

All cart models are delivered ready-assembled.

Technical details | Dimensions

Carrier profile height	780/900mm
Dimensions 	External dimensions of cart with castors Ø 75 mm (WxDxH): <ul style="list-style-type: none"> · 575x575x982 mm · 575x575x1102 mm
	External dimensions shelves (metal): <ul style="list-style-type: none"> · 500x390 mm (WxD) · 280x400 mm (WxD) · 319/290x225 mm (W1/W2xD)
	External dimensions storage shelf (metal) (W1/W2xD): 312/284x400 mm
	Drawers (metal) (WxDxH): External dimensions: 373x260x130 mm Inner dimensions 310x230x109 mm
Loading capacity	Total loading capacity: up to approx. 25 kg
	Weight-bearing capacity shelves (metal): up to approx. 5 kg/20 kg
	Weight-bearing capacity drawers (metal): up to approx. 5 kg

Materials | Colours

Carrier profile	Aluminium, anodised
Profile cover	Profile cap, stainless steel, polished
Chassis/cart	Cast aluminium, finely structured, matte, powder-coated, basalt grey (RAL 7012)
Decorative stripes on profile	<p>PVC:</p> <ul style="list-style-type: none"> · grey (6502B as per NCS) · blue (3050-R80B as per NCS) <p>Other colour combinations or customer-specific colours available on request.</p>
Shelves	Steel, powder-coated, white aluminium (RAL 9006) or ABS plastic, orion grey
Shelf with tray	Plastic tray (polycarbonate), white, metal shelf
Push-handle	Steel, powder-coated, basalt grey (RAL 7012) (2 side push-handles on the shelf)
Drawer handle	<p>Plastic design handles, grey as standard, 3 trend colours:</p> <ul style="list-style-type: none"> · grey (6502B as per NCS) · blue (3050-R80B as per NCS) · yellow (0050-Y as per NCS) <p>Other colour combinations or customer-specific colours available on request.</p>
Drawer (metal)	Steel, powder-coated, white aluminium (RAL 9006)
Castors	<p>High-quality, easy-swivel castors Ø front/back:</p> <ul style="list-style-type: none"> · polyamide housing, non-marking, grey polyurethane tyres · 4 castors with wall protection and 2 with locking device · 75 mm as standard · 100mm optional <p>Electrically conductive castors: optionally available for a surcharge</p>

Configuration examples



Ref. No. 27011

bravo® basic cart 1

Configuration	<ul style="list-style-type: none"> · carrier profile 780mm in height · 1 metal shelf 2, 280x400 mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x985
Ref. No. Decorative stripes – grey – blue	27011 27012



Ref. No. 27021

bravo® basic cart 2

Configuration	<ul style="list-style-type: none"> · carrier profile 780mm in height · 1 metal shelf 3, 444x400 mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 2 side push-handles, basalt grey (RAL 7012) · 1 chassis with 4 castors Ø 75 mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x985
Ref. No. Decorative stripes – grey – blue	27021 27022



Ref. No. 27031

bravo® basic cart 3

Configuration	<ul style="list-style-type: none"> · carrier profile 780mm in height · 1 metal shelf top, 280x400 mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 1 metal shelf 4 below, shelf 319/290x225 (W1/B2xD) · 1 chassis with 4 castors Ø 75 mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x985
Ref. No. Decorative stripes – grey – blue	27031 27032

bravo® basic cart 5


Ref. No. 27071

Configuration	<ul style="list-style-type: none"> · carrier profile 780mm in height · 1 metal shelf top, 280x400 mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 1 shelf with tray: <ul style="list-style-type: none"> – metal shelf, 324x240mm (WxD) – polycarbonate tray for accessories, pull-out, white · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x985
Ref. No. Decorative stripes – grey – blue	27071 27072

bravo® basic cart 6


Ref. No. 27361

Configuration	<ul style="list-style-type: none"> · carrier profile 780mm in height · 1 metal shelf 1, 500x390mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x985
Ref. No. Decorative stripes – grey – blue	27361 27362

bravo® basic cart 8


Ref. No. 27561

Configuration	<ul style="list-style-type: none"> · carrier profile 780mm in height · 1 metal shelf top, 500x390mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 1 shelf with metal drawer with handle · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x985
Ref. No. Decorative stripes – grey – blue	27561 27562

Configuration examples



Ref. No. 27531

bravo® basic cart 9

Configuration	<ul style="list-style-type: none"> · carrier profile 780mm in height · 1 shelf with metal drawer with handle · 2 stainless steel profile covers · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x985
Ref. No. Decorative stripes – grey – blue	27531 27532



Ref. No. 27501

bravo® equipment cart 1

Configuration	<ul style="list-style-type: none"> · carrier profile 900mm in height · 1 metal shelf 1 top, 500x390mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 1 shelf with tray: <ul style="list-style-type: none"> – metal shelf, 324x240mm (WxD) – polycarbonate tray for accessories, pull-out, white · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x1105
Ref. No. Decorative stripes – grey – blue	27501 27502



Ref. No. 27191

bravo® equipment cart 2

Configuration	<ul style="list-style-type: none"> · carrier profile 900mm in height · 1 metal shelf 1 top, 500x390mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 1 storage basket, to front of carrier profile · 1 outlet strip (3-socket), front between the carrier profiles · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x1105
Ref. No. Decorative stripes – grey – blue	27191 27192



Ref. No. 27151

bravo® computer cart 1

Configuration	<ul style="list-style-type: none"> · carrier profile 780mm in height · 1 metal shelf 2, 280x400mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 1 TFT mount, VESA standard 75/100mm, not height-adjustable, rotatable and tiltable · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x1392 (with TFT mount)
Ref. No. Decorative stripes – grey – blue	27151 27152



Ref. No. 27571

bravo® tablet cart 3

Configuration	<ul style="list-style-type: none"> · carrier profile 780mm in height · 1 metal shelf 1, 500x390mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 1 tablet holder with flex-shaft, secure attachment of a tablet/iPad (10"–15"), max. tablet thickness 12mm, load capacity up to approx. 0.8kg · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x approx. 1330 with tablet holder (without 985)
Ref. No. Decorative stripes – grey – blue	27571 27572



Ref. No. 27581

bravo® tablet cart 4

Configuration	<ul style="list-style-type: none"> · carrier profile 900mm in height · 1 tablet holder with flex-shaft, secure attachment of a tablet/iPad (10"–15"), max. tablet thickness 12mm, load capacity up to approx. 0.8kg · all-around quad push-handle, powder-coated aluminium (RAL 9006), can also be used as device rail 10x25 (DxH) for accessories, load capacity up to approx. 5kg · 1 notch at the back for cabling · 1 storage basket, to front of carrier profile · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x approx. 1450 with tablet holder (without 1105)
Ref. No. Decorative stripes – grey – blue	27581 27582

bravo®

Configuration examples



Ref. No. 27181

bravo® equipment cart ER 1

Configuration	<ul style="list-style-type: none"> · carrier profile 780mm in height · 1 metal shelf 2, 280x400mm (WxD) · 1 pull-grip · 1 notch at the back for cabling · 1 aluminium cable and hose holder, with push-button adjustment, radial range 950mm · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575xapprox. 985 (without cable and hose holder)
Ref. No. Decorative stripes – grey – blue	27181 27182



Ref. No. 27471

bravo® hygiene cart 1

Configuration	<ul style="list-style-type: none"> · carrier profile 900mm in height · 1 disinfectant dispenser for 1-litre bottles · 1 drip tray, removable · 1 chassis with 4 castors Ø 75mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x1105
Ref. No. Decorative stripes – grey – blue	27471 27472



Ref. No. 27081

bravo® hygiene cart 2

Configuration	<ul style="list-style-type: none"> · carrier profile 900 mm in height · 1 dispensing box holder (single) for disposal gloves or tissues · 1 disinfectant dispenser for 1-litre bottles · 1 drip tray, removable · 1 chassis with 4 castors Ø 75 mm, 2 lockable, with wall protectors
Dimensions in mm (WxDxH)	575x575x1246
Ref. No. Decorative stripes – grey – blue	 27081 27082



Ref. No. 27051

Chassis

Configuration	<ul style="list-style-type: none"> · 1 two-part carrier profile, aluminium · 1 cast-aluminium frame, powder-coated, basalt grey (RAL 7012), finely structured, matte · 4 castors Ø 75 mm, 2 lockable, with wall protectors · load capacity up to approx. 25 kg 	
	· carrier profile 780mm in height	· carrier profile 900mm in height
Dimensions in mm (WxDxH)	575x575x982	575x575x1102
Ref. No. Decorative stripes – grey – blue	27051 27052	27061 27062



Ref. No. 27048

Shelf 1

Configuration	<ul style="list-style-type: none"> · 1 metal shelf, powder-coated, white aluminium (RAL 9006) · 1 pull-grip · 1 notch at the back for cabling · load capacity up to approx. 20 kg 	
Dimensions in mm (WxD)	500x390	
Ref. No.	27048	



Ref. No. 27311

Shelf 2

Configuration	<ul style="list-style-type: none"> · 1 metal shelf, powder-coated, white aluminium (RAL 9006) · 1 pull-grip · 1 notch at the back for cabling · weight-bearing capacity up to approx. 20 kg 	
Dimensions in mm (WxD)	280x400	
Ref. No.	27311	



Ref. No. 27351

Shelf 3

Configuration	<ul style="list-style-type: none"> · 1 metal shelf, powder-coated, white aluminium (RAL 9006) · 1 pull-grip · 1 notch at the back for cabling · weight-bearing capacity up to approx. 20 kg 	
Dimensions in mm	External dimensions: 444x400 (WxD) Usable surface area: 312/284x400 (W1/W2xD)	
Ref. No.	27351	



Ref. No. 27038

Shelf with tray

Configuration	<ul style="list-style-type: none"> · 1 metal shelf, powder-coated, white aluminium (RAL 9006) · 1 polycarbonate, pull-out drawer for accessories, white · disinfectant-safe · dishwasher safe · temperature resistant from -30° to +120 °C · load capacity of shelf up to approx. 3 kg/tray up to approx. 2 kg
Dimensions in mm (WxDxH)	Dimensions shelf: 324x240x76 Inner dimensions tray: 230x150x55
Ref. No.	27038



Ref. No. 27241

Shelf with drawer

Configuration	<ul style="list-style-type: none"> · 1 metal shelf, 373x260 mm (WxD), powder-coated, white aluminium (RAL 9006) · 1 metal drawer · 1 handle · load capacity from the shelf up to approx. 5 kg · load capacity from the drawer up to approx. 5 kg
Dimensions in mm (WxDxH)	External dimensions: 373x260x130 Inner dimensions: 310x230x109
Ref. No. Drawer handles – grey – blue – yellow	27241 27242 27244



Ref. No. 27251

Under-mounting drawer

Configuration	<ul style="list-style-type: none"> · 1 metal drawer, powder-coated, white aluminium (RAL 9006) · 1 handle · load capacity up to approx. 5 kg Note: Drawer without shelf.
Inner dimensions in mm (WxDxH)	310x230x109
Ref. No. Drawer handles – grey – blue – yellow	27251 27252 27254

System components



Ref. No. 27391

Shelf 4

Configuration

- 1 metal shelf, powder-coated, white aluminium (RAL 9006)
- weight-bearing capacity up to approx. 5 kg

For mounting on carrier profile.

Dimensions

in mm (WxD)

Usable surface area: 280x200
 External dimensions: 319/290x225 (W1/W2xD)

Ref. No.

27391



Ref. No. 27005

Shelf with cable and hose holder

Configuration

- 1 metal shelf, powder-coated, white aluminium (RAL 9006)
- 1 pull-grip
- silver-grey, aluminium tubing
- with push-button adjustment
- 2 plastic joints
- 2 cable-routing elements
- including mounting materials
- radial range approx. 950mm
- load capacity of shelf up to approx. 20kg

Dimensions

in mm (WxD)

External dimensions: 280x400
 Usable surface area: 280x328

Ref. No.

27005



Ref. No. 27004

Shelf with TFT mount

Configuration

- 1 metal shelf, powder-coated, white aluminium (RAL 9006)
- 1 pull-grip
- 1 TFT mount, VESA standard 75/100 mm, on a rod, not height-adjustable, rotatable and tiltable
- for attaching flat monitors
- median height 348 mm to shelf
- load capacity from the shelf up to approx. 20kg
- load capacity from the TFT mount up to approx. 10kg

Dimensions

in mm (WxD)

External dimensions: 280x400
 Usable surface area: 280x325

Ref. No.

27004



Ref. No. 27986

Castors, Ø 75 mm

Configuration	<ul style="list-style-type: none"> · quality, easy-swivel castors, double-castors · additional charge for 4 castors Ø 75 mm · electrically conductive · with wall protectors · 2 castors with locking mechanism
Ref. No.	27986



Ref. No. 27907

Castors, Ø 100 mm

Configuration	<ul style="list-style-type: none"> · quality, easy-swivel castors, double-castors · additional charge for 4 castors Ø 100 mm · with wall protectors · 2 castors with locking mechanism 	
	—	· electrically conductible
Ref. No.	27907	27906



Ref. No. 19976

Device rail front

Configuration	<ul style="list-style-type: none"> · stainless steel, polished · for attachment to front of carrier profile · with two mounting points · load capacity 10kg
Dimensions in mm (WxDxH)	336 x 10 x 25 Note: Dimensions without mounting material.
Ref. No.	19976



Ref. No. 27046

Storage basket

Configuration	<ul style="list-style-type: none"> · steel, powder-coated · white aluminium (RAL 9006) · for attachment to front of back of carrier profile · mesh width 20 mm · load capacity up to approx. 5 kg
Dimensions in mm (WxDxH)	377 x 196 x 115 (With spacer, 214 mm deep)
Ref. No.	27046