

fuego



WAARDENBURG
MEDICAL



fuego

The flexible cart for small device combinations.

Flexible functionality and outstanding design make the fuego stand out from the crowd. Since it is very manoeuvrable and sturdy, fuego is perfect for mobile application with small device combinations. It offers a large storage space over a small surface area. With its sleek design, fuego is the perfect match for a multitude of areas. fuego provides mobile solutions in offices or conference rooms, too. The fuego equipment cart series stands for dependable brand quality. It is protected with a European design patent.

Are you looking for a cart that's the perfect fit? On the next pages, we'll show you a range of particularly popular configurations for our fuego. No more and no less. The fuego truly is a unique piece that can be put together and equipped completely according to your needs.



Learn about the fuego in more detail.

Features

- Single-column equipment cart system, light-duty.
- Robust chassis with a single-column aluminium profile for attaching outlet strips.
- Integrated cable duct.
- With ergonomic push-handle as standard and cable stowing clip.
- Lifting function with freely adjustable height for more working comfort.
- Optional design in various decorative colours.

Potential function

Equipment cart | ECG cart | beauty/cosmetic cart | laptop or computer cart | ergometry cart | infusion cart | Lift cart for work while sitting or standing, with free height adjustment ...

Potential utilisation

Doctors' surgeries | clinics | cosmetic studios | offices | consultation and conference rooms ...

The fuego system components

- Robust, very manoeuvrable chassis with high-quality castors; front castors lockable with wall protectors.
- High-quality plastic shelves with moulded edge and recessed grip or metal shelves.
- Custom dimensions on request.
- Elegant drawer made of metal or plastic, easy-glide telescopic rails on ball-bearings (several drawers under one shelf possible).
- With large ergonomic push-handle as standard, high-quality powder coating in white aluminium (RAL 9006).
- 2 decorative colours, custom colours for a surcharge.
- Balancing weight plate depending on the configuration of the equipment cart system to increase stability optionally available (Ref. No. 19369).


POAG equipotential bonding system conforming to DIN 42801

For your safety, selected equipment carts are equipped ex works with a POAG pin. The POAG pin (Ref. No. 19529) is attached to the outside of the rear of the carrier profile. Retrofitting by qualified specialist personnel on the part of the customer is also possible.

Note

All cart models are delivered ready-assembled.

Technical details | Dimensions

Carrier profile height	800/950mm
Dimensions	<p>External dimensions of cart with castors Ø 75/100 mm (WxDxH):</p> <ul style="list-style-type: none"> · 545x572x987 mm · 545x572x1137 mm
	<p>ABS shelves: External dimensions: 500x390x53 mm (WxDxH) Usable surface area: 476/432x340 mm (W1/W2xD)</p> <p>Metal shelves: External dimensions: 495x387 mm (WxD) Usable surface area: 495/448x369/330 mm (W1/W2xD1/D2)</p>
Loading capacity	<p>Total loading capacity: up to approx. 60kg</p> <hr/> <p>Weight-bearing capacity ABS shelves: up to approx. 15 kg/25kg</p> <p>Weight-bearing capacity metal shelves: up to approx. 25kg</p> <hr/> <p>Weight-bearing capacity ABS drawers: up to approx. 5kg</p> <p>Weight-bearing capacity metal drawers: up to approx. 5kg</p>

Materials | Colours

Carrier profile	Aluminium, anodised
Profile cover	Plastic, ABS, grey (6502B as per NCS)
Chassis/cart	Cast aluminium, finely structured, matte, powder-coated, basalt grey (RAL 7012)
Decorative stripes on profile	<p>PVC:</p> <ul style="list-style-type: none"> · grey (6502B as per NCS) · blue (3050-R80B as per NCS) <p>Other colour combinations or customer-specific colours available on request.</p>
Shelves	Steel, powder-coated, white aluminium (RAL 9006) or ABS plastic, orion grey
Rail for plates	Tubular steel, powder-coated, white aluminium (RAL 9006)
Push-handle	Tubular steel, powder-coated, white aluminium (RAL 9006)
Drawer handle	<p>Plastic design handles, grey as standard, 3 trend colours:</p> <ul style="list-style-type: none"> · grey (6502B as per NCS) · blue (3050-R80B as per NCS) · yellow (0050-Y as per NCS) <p>Other colour combinations or customer-specific colours available on request.</p>
Drawer	ABS plastic orion grey, surface polished or metal, powder-coated, white aluminium (RAL 9006)
Castors	<p>High-quality, easy-swivel castors Ø front/back:</p> <ul style="list-style-type: none"> · polyamide housing, non-marking, grey polyurethane tyres · 75 mm/ 100 mm as standard · 100 mm/ 125 mm optional <p>Electrically conductive castors: optionally available for a surcharge</p>



Ref. No. 30241

fuego basic cart 1, metal

Configuration	<ul style="list-style-type: none"> · carrier profile 800mm in height · 1 metal shelf, 495x387 mm (WxD) · powder-coated · white aluminium (RAL 9006) · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75 mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x987
Ref. No. Decorative stripes – grey – blue	30241 30242



Ref. No. 30001

fuego basic cart 1, ABS

Configuration	<ul style="list-style-type: none"> · carrier profile 800mm in height · 1 ABS plastic shelf, 500x390 mm (WxD) · with moulded edge and recessed grip · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75 mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x987
Ref. No. Decorative stripes – grey – blue	30001 30002



Ref. No. 30231

fuego basic cart 2, metal

Configuration	<ul style="list-style-type: none"> · carrier profile 800mm in height. · 2 metal shelves, 495x387 mm (WxD) · powder-coated · white aluminium (RAL 9006), · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75 mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x987
Ref. No. Decorative stripes – grey – blue	30231 30232



Ref. No. 30031

fuego basic cart 2, ABS

Configuration	<ul style="list-style-type: none"> · carrier profile 800mm in height · 2 ABS plastic shelves, 500x390 mm (WxD) · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75 mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x987
Ref. No. Decorative stripes – grey – blue	30031 30032

fuego

Configuration examples



Ref. No. 30041

fuego basic cart 3, ABS

Configuration	<ul style="list-style-type: none"> · carrier profile 950 mm in height. · 3 ABS plastic shelves, 500x390 mm (WxD) · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100 mm back, Ø 75 mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x1137
Ref. No. Decorative stripes – grey – blue	30041 30042



Ref. No. 30221

fuego basic cart 1 with drawer, metal

Configuration	<ul style="list-style-type: none"> · carrier profile 800 mm in height · 1 metal shelf, 500x400 mm (WxD), usable storage area 450x380 mm (WxD) · 1 metal drawer · powder-coated · white aluminium (RAL 9006) · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100 mm back, Ø 75 mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x987
Ref. No. Decorative stripes – grey – blue	30221 30222



Ref. No. 30011

fuego basic cart 1 with drawer, ABS

Configuration	<ul style="list-style-type: none"> · carrier profile 800mm in height · 1 ABS plastic shelf, 500x390mm (WxD) · 1 drawer ABS plastic · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x987
Ref. No. Decorative stripes – grey – blue	30011 30012



Ref. No. 30021

fuego basic cart 1 with storage basket, ABS

Configuration	<ul style="list-style-type: none"> · carrier profile 800mm in height · 1 ABS plastic shelf, 500x390mm (WxD) · 1 pull-out storage basket (Ref. No. 19268), 418x268x100mm (WxDxH) · steel, powder-coated · white aluminium (RAL 9006) · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x987
Ref. No. Decorative stripes – grey – blue	30021 30022



Ref. No. 30051

fuego laptop cart, metal

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 950mm in height · 1 laptop shelf, usable area, 491 x 370mm (WxD) · 1 mouse tray, pull-out from left or right, 130x160mm (WxD) · 1 rear storage basket (Ref. No. 19146), 377 x 196 x 115 mm (WxDxH), load capacity 5 kg (evenly distributed weight) · steel, powder-coated · white aluminium (RAL 9006) · 1 outlet strip (4-socket), front on carrier profile · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x1137
Ref. No. Decorative stripes – grey – blue	30051 30052



Ref. No. 30341

fuego tablet cart, ABS

Configuration	<ul style="list-style-type: none"> · 1 carrier profile 800mm in height · 1 ABS plastic shelf, 500 x 390mm (WxD) · 1 tablet holder with flex-shaft, secure attachment of a tablet/iPad (10"-15"), max. tablet thickness 12 mm, load capacity up to approx. 0.8kg · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x1320
Ref. No. Decorative stripes – grey – blue	30341 30342



Ref. No. 30311

fuego equipment cart MP 1, metal

Configuration	<ul style="list-style-type: none"> · 1 carrier profile 800mm in height · 2 metal shelves with handle, 350x400mm (WxD) · steel, powder-coated · white aluminium (RAL 9006) · 1 outlet strip (3-socket), front on carrier profile · 1 TFT mount, VESA standard 75/100mm, height-adjustable, tiltable · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x1127
Ref. No. Decorative stripes – grey – blue	30311 30312



Ref. No. 19751

fuego computer cart 1, ABS

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 800mm in height · 1 ABS plastic shelf, 500x390mm (WxD) · 1 outlet strip (4-socket), front on carrier profile · 1 alegro monitor bracket, VESA standard 75/100mm, height-adjustable, rotatable and tiltable · 1 pull-out keyboard/mouse tray · steel, powder-coated · white aluminium (RAL 9006) · 1 PC mount, right, load capacity up to approx. 8kg · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable · 2 wall protectors <p>PC mount can be assembled on the left, if desired. Please specify when ordering.</p>
Dimensions in mm (WxDxH)	545x572x1400
Ref. No. Decorative stripes – grey – blue	19751 19752



Ref. No. 19581

fuego computer cart 2, ABS

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 800mm in height · 1 ABS plastic shelf, 500x390mm (WxD) · 1 outlet strip (4-socket), front on carrier profile · 1 TFT mount, VESA standard 75/100mm, height adjustable, tiltable · 1 pull-out keyboard/mouse tray · steel, powder-coated · white aluminium (RAL 9006) · 1 PC mount, right, weight-bearing capacity up to approx. 8kg · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable · 2 wall protectors <p>PC mount can be assembled on the left, if desired. Please specify when ordering.</p>
Dimensions in mm (WxDxH)	545x572x1457 (incl. TFT mount)
Ref. No. Decorative stripes – grey – blue	19581 19582



Ref. No. 19621

fuego equipment cart ER 1, ABS

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 950mm in height · 2 ABS plastic shelves, 500x290mm (WxD) · 1 front storage basket (Ref. No. 19135), 418/374x268x100mm (B1/B2xDxH) · steel, powder-coated · white aluminium (RAL 9006) · 1 aluminium cable and hose holder, with push-button adjustment, radial range approx. 950mm · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	575x572x1137 (without cable and hose holder)
Ref. No. Decorative stripes – grey – blue	19621 19622



Ref. No. 30171

fuego equipment cart EC 1, ABS

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 800mm in height · 1 mount for a Strässle suction system DT 80/100 · 2 ABS plastic shelves, 500x390 mm (WxD) · 1 TFT mount with sockets, VESA standard 75/100mm, height-adjustable, not rotatable or tiltable · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75 mm front, 2 lockable · 2 wall protectors · 1 counterweight plate, 10.5 kg Delivery WITHOUT suction system.
----------------------	--

Dimensions in mm (WxDxH)	545x572x1238
------------------------------------	--------------

Ref. No. Decorative stripes – grey – blue	30171 30172
---	----------------



Ref. No. 30361

fuego equipment cart EC 1.1, ABS

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 800mm in height · 1 mount for a for Strässle suction system · 2 ABS plastic shelves, 500x390 mm (WxD) · 1 outlet strip (4-socket), front on carrier profile · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75 mm front, 2 lockable · 2 wall protectors · 1 counterweight plate, 10.5 kg Delivery WITHOUT suction system.
----------------------	---

Dimensions in mm (WxDxH)	545x572x1051
------------------------------------	--------------

Ref. No. Decorative stripes – grey – blue	30361 30362
---	----------------



Ref. No. 30181

fuego equipment cart EC 3, ABS

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 800mm in height · 1 mount for a KISS suction system · 2 ABS plastic shelves, 500x290 mm (WxD) · 1 TFT mount, VESA standard 75/100 mm, rotatable and tiltable, height-adjustable up to 395 mm max. · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100 mm back, Ø 75 mm front, 2 lockable · 2 wall protectors · 1 counterweight plate, 10.5 kg <p>Delivery WITHOUT suction system.</p>
Dimensions in mm (WxDxH)	545x572x1524 (incl. TFT mount, not incl. suction system)
Ref. No. Decorative stripes – grey – blue	30181 30182



Ref. No. 30081

fuego equipment cart EN 1, ABS

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 950mm in height · 1 Endoscope holder, rod made of stainless steel, holding plate made of plexiglass, height adjustable from 1349 to max. 2199 mm · 3 ABS plastic shelves, 500x290 mm (WxD) · 1 outlet strip (4-socket), front on carrier profile · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100 mm back, Ø 75 mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x1349 (with endoscope holder)
Ref. No. Decorative stripes – grey – blue	30081 30082

fuego special-purpose cart IF 1


Ref. No. 30191

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 800mm in height · 1 stainless steel infusion bottle holder, height-adjustable from min. 180 to max. 1030mm, load capacity up to approx. 4kg · 1 front storage basket (Ref. No. 19135), 418/374x268x100mm (B1/B2xDxH) · steel, powder-coated · white aluminium (RAL 9006) · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable · 2 wall protectors
Dimensions in mm (WxDxH)	545x572x1365 (with infusion bottle holder)
Ref. No. Decorative stripes – grey – blue	30191 30192

fuego lift cart 1 laptop, ABS


Ref. No. 19781

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 770mm in height · 1 ABS plastic shelf, 500x290mm (WxD) · 1 front handle with integrated release lever for freely adjustable height from 700–1000mm, for working while seated or standing · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable · 2 wall protectors · the lift unit has a load capacity of up to approx. 10kg <p>Top-mounted isolating transformer available on request. It is not practicable to mount an isolating transformer under the chassis.</p>
Dimensions in mm (WxDxH)	545x572x956
Ref. No. Decorative stripes – grey – blue	19781 19782



Ref. No. 30211

fuego lift cart 2 laptop, metal

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 770mm in height · 1 pull-out keyboard and mouse tray, usable surface area (WxD): approx. 456x325mm · 1 front handle with integrated release lever for freely adjustable height from 700–1000mm · 1 metal printer shelf, 400x370mm (WxD) · steel, powder-coated, white aluminium (RAL 9006) · 1 outlet strip (4-socket), on rear of carrier profile · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable, 2 wall protectors · the lift unit has a weight-bearing capacity of up to approx. 6kg <p>Top-mounted isolating transformer available on request. It is not practicable to mount an isolating transformer under the chassis.</p>
Dimensions in mm (WxDxH)	545x572x960
Ref. No. Decorative stripes – grey – blue	30211 30212
Surcharge for lifting unit up to approx. 15kg load capacity	19757



Ref. No. 30351

fuego lift cart 3 laptop, metal

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 770mm in height · 1 keyboard and mouse tray, usable surface area (WxD): approx. 456x325mm · 1 front handle with integrated release lever for freely adjustable height from 700–1000mm · steel, powder-coated, white aluminium (RAL 9006) · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable, 2 wall protectors · the lift unit has a load capacity of up to approx. 6kg <p>Top-mounted isolating transformer available on request. It is not practicable to mount an isolating transformer under the chassis.</p>
Dimensions in mm (WxDxH)	545x572x960
Ref. No. Decorative stripes – grey – blue	30351 30352
Surcharge for lifting unit up to approx. 15kg weight-bearing capacity	19757

fuego lift cart PC, metal


Ref. No. 19791

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 770mm in height · 1 keyboard and mouse tray, usable surface area (WxD): approx. 456x325mm · steel, powder-coated, white aluminium (RAL 9006) · 1 mouse holder, right or left · 1 front handle with integrated release lever for freely adjustable height from 700–1000mm · 1 alegro monitor bracket, load capacity up to approx. 12kg · 1 PC mount, right, load capacity up to approx. 8kg · 1 outlet strip (4-socket), on rear of carrier profile · 1 cable clip · 1 push-handle · 1 chassis with four castors, Ø 100mm back, Ø 75mm front, 2 lockable, 2 wall protectors · the lift unit has a load capacity of up to approx. 4kg <p>Top-mounted isolating transformer available on request. It is not practicable to mount an isolating transformer under the chassis. PC mount can be assembled on the left, if desired. Please specify when ordering.</p>
----------------------	---

Dimensions in mm (WxDxH)	545x572x1262
------------------------------------	--------------

Ref. No. Decorative stripes – grey – blue	19791 19792
Surcharge for lifting unit up to approx. 15 kg weight-bearing capacity	19757

fuego lift cart MPC, metal


Ref. No. 30201

Configuration	<ul style="list-style-type: none"> · 1 POAG pin · 1 carrier profile 770mm in height · 1 keyboard and mouse tray, usable surface area (WxD): approx. 456x325mm · steel, powder-coated, white aluminium (RAL 9006) · 1 mouse holder, right or left · 1 front handle with integrated release lever for freely adjustable height from 740–1040mm · 1 alegro monitor bracket, load capacity up to approx. 12kg · 1 front storage basket (Ref. No. 19135), 418/374 x 268 x 100 (B1 / B2 x D x H) · 1 cable clip · 1 push-handle · 1 chassis with four large, electrically conductive castors, Ø 125mm back, Ø 100mm front, 2 lockable, 2 wall protectors · the lift unit has a load capacity of up to approx. 9kg <p>Top-mounted isolating transformer available on request. It is not practicable to mount an isolating transformer under the chassis.</p>
----------------------	---

Dimensions in mm (WxDxH)	545x572x1305
------------------------------------	--------------

Ref. No. Decorative stripes – grey – blue	30201 30202
---	----------------



Ref. No. 19651

Chassis

Configuration	<ul style="list-style-type: none"> · 1 chassis including carrier profile · 1 push-handle · 1 cable clip · 4 castors, Ø 100 mm back, Ø 75 mm front · 2 front castors lockable with wall protectors · load capacity up to approx. 60 kg 	
	POAG Pin optionally available, see electrical accessories.	
	· 1 carrier profile 800 mm in height	· 1 carrier profile 950 mm in height
Dimensions in mm (WxDxH)	545x572x987	545x572x1137
Ref. No. Decorative stripes – grey – blue	19651 19652	19661 19662



Ref. No. 19865

Extra-long push-handle

Configuration	<ul style="list-style-type: none"> · 1 push-handle, extra long · made of tubular steel · powder-coated · white aluminium (RAL 9006) 	
Dimensions in mm (WxDxH)	445x145x276	
Ref. No.	19865	

Castors, Ø 75 mm/Ø 100 mm



Ref. No. 19481

Configuration	<ul style="list-style-type: none"> · quality, easy-swivel castors, double-castors · surcharge for 4 electrically conductive castors · 2 castors, Ø 75 mm, lockable · 2 back castors Ø 100 mm 	
Ref. No.	19481	

Castors, Ø 100 mm/Ø 125 mm



Ref. No. 19850

Configuration	<ul style="list-style-type: none"> · quality, easy-swivel castors · surcharge for 4 large, electrically conductive castors · 2 castors, Ø 100 mm, lockable · 2 back castors Ø 125 mm 	
Ref. No.	19850	



Ref. No. 19033

Shelf, ABS plastic


Configuration	<ul style="list-style-type: none"> · ABS plastic shelf, orion grey · moulded edge and recessed grip 	
	<ul style="list-style-type: none"> · load capacity up to approx. 15kg 	<ul style="list-style-type: none"> · load capacity up to approx. 25kg
Dimensions in mm	Usable surface area (W1 /W2 xD) 476/432x340 External dimensions (WxDxH): 500x390x53	
Ref. No.	19033	19237



Ref. No. 19312

Spacer

Configuration	<ul style="list-style-type: none"> · 1 spacer for ABS shelf to extend usable depth by 38mm
Ref. No.	19312



Ref. No. 19367

Front handle

Configuration	<ul style="list-style-type: none"> · 1 front handle for attaching to ABS shelves · finely structured · white aluminium (RAL 9006)
Ref. No.	19367



Ref. No. 19017

Rail

Configuration	<ul style="list-style-type: none"> · 1 rail for ABS or metal shelves · additional stability for devices · made of tubular steel · powder-coated · white aluminium (RAL 9006) · Ø 8 mm, surrounds shelf on three sides <p>Attached to profile.</p>
Ref. No.	19017



Ref. No. 19330

Shelf, metal

Configuration	<ul style="list-style-type: none"> · 1 metal shelf · powder-coated · white aluminium (RAL 9006) · fully assembled · load capacity up to approx. 25 kg <p>Custom dimensions on request.</p>
Dimensions in mm	Usable surface area (W1 /W2xD1 /D2) 495/448/369/330 External dimensions (WxD): 495x387
Ref. No.	19330



Ref. No. 19539

Shelf with handle, metal

Configuration	<ul style="list-style-type: none"> · 1 metal shelf · powder-coated · white aluminium (RAL 9006) · 1 handle · load capacity up to approx. 25 kg
Dimensions in mm (WxD)	350x400
Ref. No.	19539



Ref. No. 19951

Printer shelf

Configuration	<ul style="list-style-type: none"> · 1 metal shelf for printer · powder-coated · white aluminium (RAL 9006) · load capacity up to approx. 25 kg <p>The shelf can be attached to the side or in the middle of the carrier profile.</p>
Dimensions in mm (WxD)	430x445
Ref. No.	19951



Ref. No. 19995

Special shelf

Configuration	<ul style="list-style-type: none"> · 1 metal shelf, can be mounted recessed between the carrier profile and the chassis · powder-coated · white aluminium (RAL 9006) · load capacity up to approx. 10 kg
Dimensions in mm (WxD)	418x370
Ref. No.	19995

Laptop shelf



Ref. No. 19320

Configuration	<ul style="list-style-type: none"> · 1 metal shelf for laptops · 1 metal mouse tray, extendible to 130 mm to left or right · powder-coated · white aluminium (RAL 9006) · load capacity up to approx. 25 kg
Dimensions in mm	Dimensions of usable area (WxD): 491x370 Mouse tray (WxD): 130x160
Ref. No.	19320

Drawer, ABS plastic



Ref. No. 19300



Ref. No. 19376

Configuration	<ul style="list-style-type: none"> · 1 drawer, ABS plastic · orion grey, polished surface · easy-glide, ball-bearing mounted telescopic guides · click-in mechanism · load capacity up to approx. 5 kg <p>For fixing beneath an ABS shelf or under an ABS drawer.</p>
Dimensions in mm (WxDxH)	External dimensions: 494 x 386 x 131 Inner dimensions 409 x 314 x 94
Ref. No.	19300
Lock separate for integration into a drawer	19376

Drawer insert



Ref. No. 19597

Configuration	<ul style="list-style-type: none"> · 1 drawer insert with variable partitions · with 1 back-front and 4 side-to-side dividers · made of transparent plastic · removable · dishwasher safe <p>Only for the ABS plastic drawer available.</p>
Ref. No.	19597

Drawer for lift carts, ABS plastic



Ref. No. 30096

Configuration	<ul style="list-style-type: none"> · 1 stationary drawer on lift cart · 1 ABS plastic shelf · 1 drawer mount at bottom above the chassis, fixed to profile, not to the liftable part · load capacity up to approx. 5 kg
Dimensions in mm	Footprint dimensions of the ABS plastic shelf (W1/W2xD): 476/432x340 External dimensions of the drawer (WxDxH): 497x502x174 Inner dimensions of the drawer (WxDxH): 409x314x94
Ref. No.	30096



Ref. No. 30097

2 drawers for lift carts, ABS plastic

Configuration	<ul style="list-style-type: none"> · 2 stationary drawers on lift cart · 2 ABS plastic shelves · 1 drawer mount at bottom above the chassis, fixed to profile, not to the liftable part · load capacity up to approx. 5 kg per drawer
Dimensions in mm	External dimensions (WxDxH): 497x502x295 Footprint of the ABS plastic shelf (W1/W2xD): 476/432x340 Inner dimensions of the drawer: 409x314x94
Ref. No.	30097



Ref. No. 30251

Drawer, metal

Configuration	<ul style="list-style-type: none"> · 1 metal drawer, load capacity up to approx. 5 kg · 1 metal shelf, load capacity up to approx. 15 kg · powder-coated · white aluminium (RAL 9006) · 1 drawer handle, grey
Dimensions in mm	External dimensions of shelves (WxD): 500x400 Usable area (WxD): 450x380 External dimensions of drawer (WxDxH): 415x368x128 Inner dimensions (WxDxH): 315x320x106
Ref. No. Drawer handles – grey – blue – yellow	30251 30252 30254



Ref. No. 30261

Under-mounting drawer, metal

Configuration	<ul style="list-style-type: none"> · 1 metal drawer · powder-coated · white aluminium (RAL 9006) · 1 drawer handle, grey · load capacity up to approx. 5 kg · cascadable For assembly underneath a metal shelf or under a metal drawer.
Dimensions in mm (WxDxH)	External dimensions: 415x368x128 Inner dimensions 315x320x106
Ref. No. Drawer handles – grey – blue – yellow	30261 30262 30264



Ref. No. 19220

Keyboard / mouse tray

Configuration	<ul style="list-style-type: none"> · 1 pull-out keyboard and mouse tray · extendible to the front by 215 mm · mouse tray extendible to 100 mm to left or right · metal · powder-coated · white aluminium (RAL 9006) · fixed to carrier profile · load capacity of keyboard tray up to approx. 3 kg
Dimensions in mm	External dimensions (WxDxH): 491 x 387 x 128 Usable surface area (W1 / W2 x D1 / D2) 491 / 477 x 199 / 283
Ref. No.	19220
Mouse holder can be assembled left or right to the keyboard and mouse tray or metal shelf	19488



Ref. No. 19222

Keyboard / mouse tray including TFT mount

Configuration	<ul style="list-style-type: none"> · 1 keyboard tray with mouse pull-out · 1 TFT mount, VESA standard 75 / 100 mm, on a rod, fixed to carrier profile, height-adjustable to 330 mm, rotatable and tiltable · mouse tray extendible to 100 mm to left or right · metal · powder-coated · white aluminium (RAL 9006) · weight-bearing capacity of keyboard tray up to approx. 3 kg · load capacity of TFT mount up to approx. 10 kg
Dimensions in mm	External dimensions (WxDxH): 491 x 387 x approx. 470 Usable surface area (W1 / W2 x D1 / D2): 491 / 477 x 199 / 283
Ref. No.	19222



Ref. No. 19418

TFT mount, fixed

Configuration	<ul style="list-style-type: none"> · 1 TFT mount, VESA standard 75 / 100 mm, on a rod, fixed to carrier profile, not height-adjustable, rotatable and tiltable · steel · powder-coated · white aluminium (RAL 9006) · load capacity up to approx. 10 kg
External dimensions in mm (WxDxH)	123 x 97 x 331
Ref. No.	19418



Ref. No. 19223

TFT mount, adjustable

Configuration	<ul style="list-style-type: none"> · 1 TFT mount, VESA standard 75/100 mm, on a rod fixed to carrier profile, height-adjustable to 330 mm, rotatable and tiltable · steel · powder-coated · white aluminium (RAL 9006) · load capacity up to approx. 10 kg
External dimensions in mm (WxDxH)	123x97x411
Ref. No.	19223



Ref. No. 30007

Mount for suction system, TFT 1

Configuration	<ul style="list-style-type: none"> · 1 TFT mount with connectors for attaching Ohlenschläger suction system · internal tube diameter Ø 28.2 mm · steel · powder-coated · white aluminium (RAL 9006) · load capacity up to approx. 10 kg <p>Counterweight plate (Ref. No. 19369) required, 10.5 kg.</p>
Dimensions TFT in mm (WxDxH)	123x136x339
Ref. No.	30007



Ref. No. 30006

Mount for suction system, 1

Configuration	<ul style="list-style-type: none"> · 1 mount with connectors for attaching Ohlenschläger suction system · internal tube diameter Ø 28.2 mm · steel · powder-coated · white aluminium (RAL 9006) <p>Counterweight plate (Ref. no. 19369) required, 10.5 kg.</p>
Dimensions in mm (WxDxH)	128x96x85
Ref. No.	30006



Ref. No. 19456

Mount for suction system, adaptation

Configuration	<ul style="list-style-type: none"> · adaptation for KISS suction system · internal tube diameter Ø 19.8 mm · aluminium · load capacity up to approx. 8.5 kg <p>Counterweight plate (Ref. No. 19369) required, 10.5 kg.</p>
Dimensions in mm (WxDxH)	123x97x219
Ref. No.	19456

Mount for suction system, TFT 2



Ref. No. 19978

Configuration	<ul style="list-style-type: none"> · 1 TFT mount with connectors for attaching Strässle suction system DT 80/100 · internal tube diameter Ø 35 mm · steel · powder-coated · white aluminium (RAL 9006) · load capacity up to approx. 10 kg Counterweight plate (Ref. No. 19369) required, 10.5 kg.
Dimensions TFT in mm (WxDxH)	123x97x110
Ref. No.	19978

Mount for suction system 2



Ref. No. 19148

Configuration	<ul style="list-style-type: none"> · 1 TFT mount cap with connectors for attaching Strässle suction system DT 80/100 · internal tube diameter Ø 35 mm · steel · powder-coated · white aluminium (RAL 9006) Counterweight plate (Ref. No. 19369) required, 10.5 kg.
Dimensions TFT in mm (WxDxH)	123x135x297
Ref. No.	19148

Mount for suction system, TFT 3



Ref. No. 19981

Configuration	<ul style="list-style-type: none"> · 1 TFT mount for attaching KISS suction system · height-adjustable · internal tube diameter Ø 19.8 mm · steel · powder-coated · white aluminium (RAL 9006) · load capacity up to approx. 8.5 kg Counterweight plate (Ref. No. 19369) required, 10.5 kg.
Dimensions TFT in mm (WxDxH)	123x97x727
Ref. No.	19981

Device rail front



Ref. No. 19976

Configuration	<ul style="list-style-type: none"> · stainless steel, polished · for attachment to front of carrier profile · with two mounting points · load capacity up to approx. 10 kg
Dimensions in mm (WxDxH)	336x10x25 Note: Dimensions without mounting material.
Ref. No.	19976



Ref. No. 30016

Device rail rear

Configuration	<ul style="list-style-type: none"> · stainless steel, brushed · for attachment to rear of carrier profile · load capacity up to approx. 10kg
Dimensions in mm (WxDxH)	496x10x25 Note: Dimensions without mounting material.
Ref. No.	30016



Ref. No. 19135

Storage basket, front

Configuration	<ul style="list-style-type: none"> · steel, powder-coated · white aluminium (RAL 9006) · for attachment to front of carrier profile · mesh width 20mm · load capacity up to approx. 5kg, with even weigh distribution
Dimensions in mm (B1/B2xDxH)	418/374x268x100 (with spacer, 284mm deep)
Ref. No.	19135



Ref. No. 19146

Storage basket, rear

Configuration	<ul style="list-style-type: none"> · steel, powder-coated · white aluminium (RAL 9006) · for attachment to front of carrier profile · mesh width approx. 20mm · load capacity up to approx. 5kg, with even weigh distribution
Dimensions in mm (WxDxH)	377 x 196 x 115
Ref. No.	19146



Ref. No. 19268

Storage basket, pull-out

Configuration	<ul style="list-style-type: none"> · steel, powder-coated · white aluminium (RAL 9006) · mesh width approx. 20mm · extendible to the front by 256mm · load capacity up to approx. 5kg, with even weigh distribution
Dimensions in mm (B1/B2xDxH)	418/374x268x100 (with mounting bracket, 368mm deep)
Ref. No.	19268

Recommended accessories

You will find a larger selection of products in the accessories chapter.



Ref. No. 21696



Ref. No. 19034



Ref. No. 19035



Ref. No. 19403

Outlet socket strip, 4-socket

Configuration	<ul style="list-style-type: none"> · 4 grounding sockets · green, 2-pin on/off-switch · standard F-type (Germany) · light-grey cable · 3m-long mains cable with moulded on earthed safety plug · voltage 230V · power max. 3300VA · high-quality anodised aluminium housing <p>Note: Mounting only in the carrier profile possible. Further outlet strips are listed in the electrical accessories programme.</p>
Dimensions in mm (WxDxH)	269x44x44 mm
Ref. No.	21696

Cable and hose holder

Configuration	<ul style="list-style-type: none"> · silver-grey, aluminium tubing · push-button adjustment · 2 plastic joints · 2 cable-routing elements · including mounting materials · radial range approx. 950mm
Ref. No.	19034

Infusion bottle holder, 25 mm

Configuration	<ul style="list-style-type: none"> · stainless steel · 4 hooks · rod Ø 25 mm · 1000 mm in length · extension 850 mm · height-adjustable 180–1030 mm · for insertion into aluminium carrier profile · load capacity up to approx. 4 kg
Dimensions	Length 1000 mm
Ref. No.	19035

Holder for oxygen bottle 5/10l

Configuration	<ul style="list-style-type: none"> · for oxygen or gas bottles à 5/10l · up to Ø 145 mm · steel · powder-coated · white aluminium (RAL 9006) · including upper retaining bracket and flexible retaining strap <p>For a 10l bottle, please use an additional retaining bracket.</p>
Ref. No. Mounting right Mounting left	19403 19402